

## Filtros de jeringa de Acetato de Celulosa SCA, CHMLAB®

Ref.: CHM.SCA | Categoría: Filtros de jeringa

### Descripción

---



CHM® SCA syringe filters are designed for the quick and efficient filtration up to 100 ml of liquid. Ready-to-use units, offer high flow rates at low inlet pressures, presented in 5 pore sizes: 0.2, 0.45, 0.8, 1.2 and 5 µm and in 2 diameters; 13 and 25 mm, to fulfil your filtration requirements for clarifying/ultra cleaning. They are supplied in sterile and non-sterile versions.

#### Features:

- Hydrophilic membrane
- Low protein binding
- High throughput
- Superior strength and stability
- Up to 100 ml of sample

#### Applications:

- HPLC sample preparation
- Biological sample preparation
- Protein and enzyme filtration
- Cell Culture
- Clarification of aqueous and alcohol solutions



## Especificaciones técnicas

Referencia	Poros (µm)	Ø(mm)	Esterilidad	Presentación	Env.
SCA020015K-S	0.2	15	Sí	Box 50	1
SCA020015Q	0.2	15	No	Box 500	1
SCA020025K-S	0.2	25	Sí	Box 50	1
SCA020025H	0.2	25	No	Box 100	1
SCA020025Q	0.2	25	No	Box 500	1
SCA045015K-S	0.45	15	Sí	Box 50	1
SCA045015Q	0.45	15	No	Box 500	1
SCA045025K-S	0.45	25	Sí	Box 50	1
SCA045025H	0.45	25	No	Box 100	1
SCA045025Q	0.45	25	No	Box 500	1
SCA080015K-S	0.8	15	Sí	Box 50	1
SCA080015Q	0.8	25	No	Box 500	1
SCA080025K-S	0.8	25	Sí	Box 50	1
SCA080025H	0.8	25	No	Box 100	1
SCA080025Q	0.8	15	No	Box 500	1
SCA120015K-S	1.2	15	Sí	Box 50	1
SCA120015Q	1.2	15	No	Box 500	1
SCA120025K-S	1.2	25	Sí	Box 50	1
SCA120025H	1.2	25	No	Box 100	1
SCA120025Q	1.2	25	No	Box 500	1
SCA500015K-S	5	15	Sí	Box 50	1
SCA500015Q	5	15	No	Box 500	1
SCA500025K-S	5	25	Sí	Box 50	1
SCA500025H	5	25	No	Box 100	1
SCA500025Q	5	25	No	Box 500	1