

Filtros de membrana de Acetato de Celulosa MCA, CHMLAB®

Ref.: CHM.MCA | Categoría: Filtros de membrana

Descripción

Cellulose Acetate membranes type MCA, are recommended for the aqueous samples, biological applications and protein filtration.

Filtration membranes are composed of pure cellulose acetate that is internally supported by an inert polyester web.

Its uniform pore size and consistent flow rates ensure reliable performance.

These membranes combine high flow rates and thermal stability with very low adsorption characteristics, making the 0.2 μ m pore size perfectly suited for use in disc filter holders to sterilize aqueous solutions, buffers and media. They are also low in extractables and can be repeatedly autoclaved.

They are supplied in four pore sizes, 0.2, 0.45, 0.65 and 0.80 μ m, and in seven different diameters, 13, 25, 47, 50, 90, 142 and 293 mm (other diameters available under request).



Features:

Hydrophilic membrane

Compatible to aqueous solutions with pH 4-8, most alcohols, hydrocarbons and oils

Low extractables: ensures tests will be clean with consistent results

Exceptional dimensional strength and low binding characteristics

High flow rate

Applications:

Aqueous solutions filtration

Protein and enzyme filtration

Biological and clinical analysis

Tissue culture media sterilization



Especificaciones técnicas

Referencia	Poro (µm)	Ø(mm)	Esterilidad	Presentación	Env.
MCA020013H	0.2	13	No	Caja 100	1
MCA020025H	0.2	25	No	Caja 100	1
MCA020047H	0.2	47	No	Caja 100	1
MCA020050H	0.2	50	No	Caja 100	1
MCA020090T	0.2	90	No	Caja de 25	1
MCA020142T	0.2	142	No	Caja de 25	1
MCA020293T	0.2	293	No	Caja de 25	1
MCA045013H	0.45	13	No	Caja 100	1
MCA045025H	0.45	25	No	Caja 100	1
MCA045047H	0.45	47	No	Caja 100	1
MCA045050H	0.45	50	No	Caja 100	1
MCA045090T	0.45	90	No	Caja de 25	1
MCA045142T	0.45	142	No	Caja de 25	1
MCA045293T	0.45	293	No	Caja de 25	1
MCA065025H	0.65	25	No	Caja 100	1
MCA065047H	0.65	47	No	Caja 100	1
MCA065050H	0.65	50	No	Caja 100	1
MCA065090T	0.65	90	No	Caja de 25	1
MCA065142T	0.65	142	No	Caja de 25	1
MCA065293T	0.65	293	No	Caja de 25	1
MCA080013H	0.8	13	No	Caja 100	1
MCA080025H	0.8	25	No	Caja 100	1
MCA080047H	0.8	47	No	Caja 100	1
MCA080050H	0.8	50	No	Caja 100	1
MCA080090T	0.8	90	No	Caja de 25	1
MCA080142T	0.8	142	No	Caja de 25	1
MCA080293T	0.8	293	No	Caja de 25	1